

Intelligent Delivery...in a Flash



**LATEX
FREE**

- **EFFECTIVE** Fully integrated Smart Power system, no dose controller
- **CONVENIENT** NO SALINE NEEDED
- **FAST DELIVERY** 40 mA*min in 1:15 hrs; 80 mA*min dosage in just 2.5 hrs
- **ACCURATE** Timer-controlled, definitive shut-off at end of treatment
- **COMFORTABLE** Flexible, low-profile, no skin irritation, no adhesive burns from prolonged wear
- **COST EFFECTIVE** No costly dose controller to prime the skin; more accurate timing for precise dosage



Wireless Interactive Treatment	ActivaTek, Inc. ActivaPatch®
Power Source	6V
Interactive Intelligent Microprocessor	•
Interactive Treatment With Shut Off	•
Skin Resistance Indicator	•
Open Circuit Indicator	•
Memorizes Treatment Time	•
Dosage Options (40/80mA*min)	•
Treatment Time (80 mA*min)	2 hrs 30 min
80 mA*min Dosage	80mA*min ±8 mA*min

ActivaPatch® Specifications

Active Area	11.4 cm ²
Dosage	40 or 80 mA*min
Skin Adhesive	Hypoallergenic
Maximum Buffer Capacity	Ag/AgCl Stabilizer
Fill Volume	2.0cc
Smart Power	6V with Intelligent Microprocessor
Packaging	6 ActivaPatch Kits, Alcohol Wipes, Patient Information Cards

Distributed Exclusively by:

**PATTERSON
MEDICAL**

**SAMMONS
PRESTON**

To Order: ActivaPatch®
Sammons Preston Item #: 0815-16-608

1-800-323-5547

www.sammonspreston.com

A Product of: ActivaTek®

See our interactive website
www.activapatch.com